



## **APARAT**

**Vaga  
Balance**

## **PROIZVOĐAČ I MODEL**

**Shimadzu AX 200**

## **VAGA**

### **Kratki opis metode**

Prije vaganja, vagu je potrebno staviti u ravnotežu (to se radi pritiskom na određenu tipku). Nakon toga se predmet stavlja na sredinu vage, pričekava se da se brojevi na zaslonu stabiliziraju i očita odvaga.

### **Short description of the method**

Before weighing, balance has to be put to steady state. After that object is put in the middle of the balance, wait for the numbers on the display to stabilize and read the result.

### **Namjena**

Služi za mjerenje mase uzoraka.

### **Purpose**

Used for measuring mass of samples.

**Tehničke značajke**

Za mase od 0,1mg do 200g

**Technical characteristics**

For mass from 0.1mg to 200g

**Tip i priprava uzorka**

Anorganski/organski; prirodni/sintetski; najčešće u krutom stanju.

**Sample type and preparation**

Inorganic/organic; naturally occurring/synthetic; usually in solid state.